

Neurotrophin-3

Catalog No: #AP78384


Package Size: #AP78384-1 50ug #AP78384-2 100ug #AP78384-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Neurotrophin-3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	139-257aa
Other Names	HDNF,Nerve growth factor 2,Neurotrophic factor,NGF-2
Accession No.	P20783Gene name:NTF3
Uniprot	P20783
GeneID	4908;
Calculated MW	13.09
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Images



Background

Seems to promote the survival of visceral and proprioceptive sensory neurons.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only