

Neurotrophin-3

Catalog No: #AP78385



Package Size: #AP78385-1 50ug #AP78385-2 100ug #AP78385-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Neurotrophin-3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rat
Immunogen Description	19-139AA
Other Names	Ntf-3,HDNF,Nerve growth factor 2,Neurotrophic factor,NGF-2
Accession No.	P18280Gene name:Ntf3
Uniprot	P18280
GeneID	81737;
Calculated MW	13.31
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Seems to promote the survival of visceral and proprioceptive sensory neurons.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only