

Neurotrophin-4

Catalog No: #AP78386



Package Size: #AP78386-1 50ug #AP78386-2 100ug #AP78386-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Neurotrophin-4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rat
Immunogen Description	80-209AA
Other Names	Nt4,Ntf5,Neurotrophin-5,Neutrophic factor 4,NT-5
Accession No.	P34131Gene name:Ntf4
Uniprot	P34131
GeneID	25730;
Calculated MW	14.3
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Could serve as a target-derived trophic factor for sensory and sympathetic neurons.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only