

P-selectin

Catalog No: #AP78393

Package Size: #AP78393-1 50ug #AP78393-2 100ug #AP78393-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	P-selectin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	532-689aa
Other Names	GMRP,GRMP,CD62 antigen-like family member P,Granule membrane protein 140,Leukocyte-endothelial cell adhesion molecule 3,Platelet activation dependent granule-external membrane protein,GMP-140,LECAM3,PADGEM,CD62P
Accession No.	P16109Gene name:SELP
Uniprot	P16109
GeneID	6403;
Calculated MW	17.38
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Ca(2+)-dependent receptor for myeloid cells that binds to carbohydrates on neutrophils and monocytes. Mediates the interaction of activated endothelial cells or platelets with leukocytes. The ligand recognized is sialyl-Lewis X. Mediates rapid rolling of leukocyte rolling over vascular surfaces during the initial steps in inflammation through interaction with PSGL1.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only