

ATF-5 Polyclonal Antibody

Catalog No: #40626

Package Size: #40626-1 50ul #40626-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ATF-5 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Rt
Specificity	ATF-5 Polyclonal Antibody detects endogenous levels of ATF-5 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human ATF-5.
Target Name	ATF-5
Other Names	ATF5; ATFX; Cyclic AMP-dependent transcription factor ATF-5; cAMP-dependent transcription factor ATF-5; Activating transcription factor 5; Transcription factor ATFx
Accession No.	Swiss-Prot: Q9Y2D1NCBI Gene ID: 22809
Uniprot	Q9Y2D1
GeneID	22809;
SDS-PAGE MW	31kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/5000.

Not yet tested in other applications.

Images

Western Blot analysis of Jurkat cells using ATF-5 Polyclonal Antibody



Note: This product is for in vitro research use only