

Tumor necrosis factor

Catalog No: #AP78404



Package Size: #AP78404-1 50ug #AP78404-2 100ug #AP78404-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tumor necrosis factor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	Recombinant protein
Other Names	TNFA, TNFSF2, Cachectin, TNF-alpha, Tumor necrosis factor ligand superfamily member 2, TNF-a
Accession No.	P01375 Gene name: TNF
Uniprot	P01375
GeneID	7124;
Tag Info	His
Formulation	50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.) Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

The TNF intracellular domain (ICD) form induces IL12 production in dendritic cells.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only