

PFKFB3 antibody

Catalog No: #22584

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


| | |
|-----------------------|--|
| Product Name | PFKFB3 antibody |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | Purified by antigen-affinity chromatography. |
| Applications | WB |
| Species Reactivity | Hu |
| Immunogen Type | Peptide |
| Immunogen Description | Synthetic peptide contain a sequence corresponding to a region within amino acids 457 and 520 of PFKFB3 |
| Target Name | PFKFB3 |
| Other Names | PFK2; IPFK2 |
| Accession No. | Swiss-Prot:Q16875Gene ID:5209 |
| Uniprot | Q16875 |
| GeneID | 5209; |
| Concentration | 1mg/ml |
| Formulation | Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative. |
| Storage | Store at -20°C for long term preservation (recommended). Store at 4°C for short term use. |

Application Details

Predicted MW: 60kd

Western blotting: 1:500-1:3000

Images



Sample (30 ug of whole cell lysate)

A: A431

B: H1299

7.5% SDS PAGE

PFKFB3 antibody diluted at 1: 1000

Note: This product is for in vitro research use only