

Carcinoembryonic antigen-related cell adhesion molecule 5

Catalog No: #AP78416

Orders: order@signalwayantibody.com

Package Size: #AP78416-1 50ug #AP78416-2 100ug #AP78416-3 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Carcinoembryonic antigen-related cell adhesion molecule 5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	35-144AA
Other Names	CEA,Carcinoembryonic antigen,Meconium antigen 100,CEA,CD66e
Accession No.	P06731Gene name:CEACAM5
Uniprot	P06731
GeneID	1048;
Calculated MW	12.1
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Cell surface glycoprotein that plays a role in cell adhesion and in intracellular signaling. Receptor for E.coli Dr adhesins.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only