

Mucin-16

Catalog No: #AP78419



Package Size: #AP78419-1 50ug #AP78419-2 100ug #AP78419-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mucin-16
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	4743-4953AA
Other Names	CA125,Ovarian cancer-related tumor marker CA125,Ovarian carcinoma antigen CA125,CA-125
Accession No.	Q8WXI7Gene name:MUC16
Uniprot	Q8WXI7
GeneID	94025;
Calculated MW	23.21
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Thought to provide a protective, lubricating barrier against particles and infectious agents at mucosal surfaces.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only