

Glutamate decarboxylase 2

Catalog No: #AP78420



Package Size: #AP78420-1 50ug #AP78420-2 100ug #AP78420-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Glutamate decarboxylase 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Mouse
Immunogen Description	1-120AA
Other Names	Gad65,65 kDa glutamic acid decarboxylase,Glutamate decarboxylase 65 kDa isoform,GAD-65
Accession No.	P48320Gene name:Gad2
Uniprot	P48320
GeneID	14417;
Calculated MW	13.2
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Catalyzes the production of GABA.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only