

Serpín B4

Catalog No: #AP78421



Package Size: #AP78421-1 50ug #AP78421-2 100ug #AP78421-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Serpín B4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	1-240AA
Other Names	PI11,SCCA2,Leupin,Peptidase inhibitor 11,Squamous cell carcinoma antigen 2,PI-11,SCCA-2
Accession No.	P48594Gene name:SERPINB4
Uniprot	P48594
GeneID	6318;
Calculated MW	26.4
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

May act as a protease inhibitor to modulate the host immune response against tumor cells.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only