

# Alpha-galactosidase A

Catalog No: #AP78429



Package Size: #AP78429-1 50ug #AP78429-2 100ug #AP78429-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Alpha-galactosidase A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	32-429AA
Other Names	Alpha-D-galactosidase A, Alpha-D-galactoside galactohydrolase, Melibiase
Accession No.	P06280 Gene name: GLA
Uniprot	P06280
GeneID	2717;
Calculated MW	43.78
Tag Info	His
Formulation	50mM NaH <sub>2</sub> PO <sub>4</sub> , 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.) Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

## Background

null

## References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only