

## Alpha-2-HS-glycoprotein

Catalog No: #AP78432



Package Size: #AP78432-1 50ug #AP78432-2 100ug #AP78432-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Alpha-2-HS-glycoprotein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	19-300AA
Other Names	PRO2743,FETUA,Alpha-2-Z-globulin,Ba-alpha-2-glycoprotein,Fetuin-A
Accession No.	P02765Gene name:AHSG
Uniprot	P02765
GeneID	197;
Calculated MW	31.02
Tag Info	His
Formulation	50mM NaH <sub>2</sub> PO <sub>4</sub> , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

## Background

Promotes endocytosis, possesses opsonic properties and influences the mineral phase of bone. Shows affinity for calcium and barium ions.

## References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only