

## Fatty acid-binding protein, adipocyte

Catalog No: #AP78435



Package Size: #AP78435-1 50ug #AP78435-2 100ug #AP78435-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Fatty acid-binding protein, adipocyte
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	2-132AA
Other Names	Adipocyte lipid-binding protein,Adipocyte-type fatty acid-binding protein,Fatty acid-binding protein 4,ALBP,A-FABP
Accession No.	P15090Gene name:FABP4
Uniprot	P15090
GeneID	2167;
Calculated MW	14.41
Tag Info	His
Formulation	50mM NaH <sub>2</sub> PO <sub>4</sub> , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

## Background

Lipid transport protein in adipocytes. Binds both long chain fatty acids and retinoic acid. Delivers long-chain fatty acids and retinoic acid to their cognate receptors in the nucleus.

## References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only