

# Corticosteroid-binding globulin

Catalog No: #AP78437



Package Size: #AP78437-1 50ug #AP78437-2 100ug #AP78437-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Corticosteroid-binding globulin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Mouse
Immunogen Description	23-397AA
Other Names	Cbg, Serpin A6, Transcortin
Accession No.	Q06770 Gene name: Serpina6
Uniprot	Q06770
GeneID	12401;
Calculated MW	41.25
Tag Info	His
Formulation	50mM NaH <sub>2</sub> PO <sub>4</sub> , 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.) Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

## Background

Major transport protein for glucocorticoids and progestins in the blood of almost all vertebrate species.

## References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only