

## Insulin-like 3

Catalog No: #AP78438

Package Size: #AP78438-1 50ug #AP78438-2 100ug #AP78438-3 1mg


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Insulin-like 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	1-131aa
Other Names	RLF,RLNL,Leydig insulin-like peptide,Relaxin-like factor,Ley-I-L
Accession No.	P51460Gene name:INSL3
Uniprot	P51460
GeneID	3640;
Calculated MW	14.41
Tag Info	His
Formulation	50mM NaH <sub>2</sub> PO <sub>4</sub> , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

## Images



## Background

Seems to play a role in testicular function. May be a trophic hormone with a role in testicular descent in fetal life. Is a ligand for LGR8 receptor.

## References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

---

Note: This product is for in vitro research use only