

N-acetylgalactosamine-6-sulfatase

Catalog No: #AP78440

Package Size: #AP78440-1 50ug #AP78440-2 100ug #AP78440-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	N-acetylgalactosamine-6-sulfatase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	27-379AA
Other Names	Chondroitinsulfatase,Galactose-6-sulfate sulfatase,N-acetylgalactosamine-6-sulfate sulfatase,Chondroitinase,GalN6S,GalNAc6S sulfatase
Accession No.	P34059Gene name:GALNS
Uniprot	P34059
GeneID	2588;
Calculated MW	38.83
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

null

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only