

# Neurosecretory protein VGF

Catalog No: #AP78441



Package Size: #AP78441-1 50ug #AP78441-2 100ug #AP78441-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Neurosecretory protein VGF
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	353-447AA
Other Names	null
Accession No.	O15240Gene name:VGF
Uniprot	O15240
GeneID	7425;
Calculated MW	10.45
Tag Info	His
Formulation	50mM NaH <sub>2</sub> PO <sub>4</sub> , 500mM NaCl Buffer with 500mM Imidazole, 10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

## Background

Antimicrobial peptide VGF[554-577]: Has bactericidal activity against *M. luteus*, and antifungal activity against *P. Pastoris*.

## References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only