

Transforming growth factor beta-2

Catalog No: #AP78446



Package Size: #AP78446-1 50ug #AP78446-2 100ug #AP78446-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Transforming growth factor beta-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	21-302AA
Other Names	BSC-1 cell growth inhibitor,Cetermin,Glioblastoma-derived T-cell suppressor factor,Polyergin,G-TSF
Accession No.	P61812Gene name:TGFB2
Uniprot	P61812
GeneID	7042;
Calculated MW	31.02
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

TGF-beta 2 has suppressive effects on interleukin-2 dependent T-cell growth.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only