

Interferon beta

Catalog No: #AP78449

Package Size: #AP78449-1 50ug #AP78449-2 100ug #AP78449-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Interferon beta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	22-187AA
Other Names	IFB,IFNB,Fibroblast interferon
Accession No.	P01574Gene name:IFNB1
Uniprot	P01574
GeneID	3456;
Calculated MW	18.26
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Has antiviral, antibacterial and anticancer activities.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only