

BAF170 Polyclonal Antibody

Catalog No: #40631

Package Size: #40631-1 50ul #40631-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	BAF170 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IF ELISA
Species Reactivity	Hu Ms
Specificity	BAF170 Polyclonal Antibody detects endogenous levels of BAF170 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human BAF170.
Target Name	BAF170
Other Names	SMARCC2; BAF170; SWI/SNF complex subunit SMARCC2; BRG1-associated factor 170; BAF170; SWI/SNF complex 170 kDa subunit; SWI/SNF-related matrix-associated actin-dependent regulator of chromatin subfamily C member 2
Accession No.	Swiss-Prot: Q8TAQ2NCBI Gene ID: 6601
Uniprot	Q8TAQ2
GeneID	6601;
SDS-PAGE MW	160kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunofluorescence: 1/200 - 1/1000.


ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of COLO205 cells using BAF170 Polyclonal Antibody



Western Blot analysis of 293T AD293 22RV1 HELA HepG2 cells using BAF170 Polyclonal Antibody

Note: This product is for in vitro research use only