

Trypsin-2

Catalog No: #AP78454


Package Size: #AP78454-1 50ug #AP78454-2 100ug #AP78454-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Trypsin-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	16-247AA
Other Names	TRY2, TRYP2, Anionic trypsinogen, Serine protease 2, Trypsin II
Accession No.	P07478 Gene name: PRSS2
Uniprot	P07478
GeneID	5645;
Calculated MW	25.52
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.) Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Images



Background

In the ileum, may be involved in defensin processing, including DEFA5.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only