

Myosin-9

Catalog No: #AP78457



Package Size: #AP78457-1 50ug #AP78457-2 100ug #AP78457-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Myosin-9
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Immunogen Description	N-terminal
Other Names	Myosin heavy chain 9, Myosin heavy chain, non-muscle IIa, Non-muscle myosin heavy chain IIa, NMMHC II-a, NMMHC-IIA
Accession No.	Q258K2 Gene name: MYH9
Uniprot	Q258K2
GeneID	481280;
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.) Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Cellular myosin that appears to play a role in cytokinesis, cell shape, and specialized functions such as secretion and capping.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only