

## Propionyl-CoA carboxylase beta chain, mitochondrial

Catalog No: #AP78462



Package Size: #AP78462-1 50ug #AP78462-2 100ug #AP78462-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Propionyl-CoA carboxylase beta chain, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	C-terminal
Other Names	Propanoyl-CoA:carbon dioxide ligase subunit beta
Accession No.	P05166Gene name:PCCB
Uniprot	P05166
GeneID	5096;
Tag Info	His
Formulation	50mM NaH <sub>2</sub> PO <sub>4</sub> , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

## Background

null

## References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only