

Beta-defensin 3

Catalog No: #AP78468

Package Size: #AP78468-1 50ug #AP78468-2 100ug #AP78468-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Beta-defensin 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Mouse
Immunogen Description	23-63AA
Other Names	Bd3,Defensin, beta 3
Accession No.	Q9WTL0Gene name:Defb3
Uniprot	Q9WTL0
GeneID	27358;
Calculated MW	4.51
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Antimicrobial activity against Gram-negative bacteria E.coli and P.aeruginosa.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only