

Pentraxin-related protein PTX3

Catalog No: #AP78469

Package Size: #AP78469-1 50ug #AP78469-2 100ug #AP78469-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Pentraxin-related protein PTX3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	18-381AA
Other Names	TNFAIP5,TSG14,Pentraxin-related protein PTX3,Tumor necrosis factor alpha-induced protein 5,Tumor necrosis factor-inducible gene 14 protein,TNF alpha-induced protein 5,TSG-14
Accession No.	P26022Gene name:PTX3
Uniprot	P26022
GeneID	5806;
Calculated MW	40.04
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Plays a role in the regulation of innate resistance to pathogens, inflammatory reactions, possibly clearance of self-components and female fertility.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only