

Spectrin alpha chain, brain

Catalog No: #AP78470



Package Size: #AP78470-1 50ug #AP78470-2 100ug #AP78470-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Spectrin alpha chain, brain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	149-253AA
Other Names	SPTA2,Alpha-II spectrin,Fodrin alpha chain,Spectrin, non-erythroid alpha chain
Accession No.	Q13813Gene name:SPTAN1
Uniprot	Q13813
GeneID	6709;
Calculated MW	11.55
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Fodrin, which seems to be involved in secretion, interacts with calmodulin in a calcium-dependent manner and is thus candidate for the calcium-dependent movement of the cytoskeleton at the membrane.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only