

Cathepsin L1

Catalog No: #AP78473

Package Size: #AP78473-1 50ug #AP78473-2 100ug #AP78473-3 1mg


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Cathepsin L1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	114-333AA
Other Names	CTSL1,Cathepsin L, Major excreted protein, MEP
Accession No.	P07711 Gene name: CTSL
Uniprot	P07711
GeneID	1514;
Calculated MW	24.2
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.) Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Images



Background

Important for the overall degradation of proteins in lysosomes.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only