

Beta-2-glycoprotein 1

Catalog No: #AP78475

Package Size: #AP78475-1 50ug #AP78475-2 100ug #AP78475-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Beta-2-glycoprotein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	20-345AA
Other Names	B2G1,APC inhibitor,Activated protein C-binding protein,Anticardiolipin cofactor,Apolipoprotein H,Beta-2-glycoprotein I,Apo-H,B2GPI
Accession No.	P02749Gene name:APOH
Uniprot	P02749
GeneID	350;
Calculated MW	35.86
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Binds to various kinds of negatively charged substances such as heparin, phospholipids, and dextran sulfate. May prevent activation of the intrinsic blood coagulation cascade by binding to phospholipids on the surface of damaged cells.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only