

Cathepsin B

Catalog No: #AP78508

Package Size: #AP78508-1 50ug #AP78508-2 100ug #AP78508-3 1mg


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Cathepsin B
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	80-203AA
Other Names	CPSB,APP secretase,Cathepsin B1,APPS
Accession No.	P07858Gene name:CTSB
Uniprot	P07858
GeneID	1508;
Calculated MW	13.64
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Images



Background

Thiol protease which is believed to participate in intracellular degradation and turnover of proteins. Has also been implicated in tumor invasion and metastasis.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only