

cAMP-specific 3',5'-cyclic phosphodiesterase 4B

Catalog No: #AP78518



Package Size: #AP78518-1 50ug #AP78518-2 100ug #AP78518-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	cAMP-specific 3',5'-cyclic phosphodiesterase 4B
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rat
Immunogen Description	1-172AA
Other Names	DPDE4
Accession No.	P14646Gene name:Pde4b
Uniprot	P14646
GeneID	24626;
Calculated MW	18.92
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Hydrolyzes the second messenger cAMP, which is a key regulator of many important physiological processes.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only