

Glycoprotein hormones alpha chain

Catalog No: #AP78521

Package Size: #AP78521-1 50ug #AP78521-2 100ug #AP78521-3 1mg


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Glycoprotein hormones alpha chain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	25-116aa
Other Names	Anterior pituitary glycoprotein hormones common subunit alpha,Choriogonadotropin alpha chain,Chorionic gonadotrophin subunit alpha,Follicle-stimulating hormone alpha chain,Follitropin alpha chain,Luteinizing hormone alpha chain,Lutropin alpha chain,Thyroid-stimulating hormone alpha chain,Thyrotropin alpha chain,CG-alpha,FSH-alpha,LSH-alpha,TSH-alpha
Accession No.	P01215Gene name:CGA
Uniprot	P01215
GeneID	1081;
Calculated MW	10.12
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Images



Background

null

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only