

## Thyrotropin subunit beta

Catalog No: #AP78526



Package Size: #AP78526-1 50ug #AP78526-2 100ug #AP78526-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Thyrotropin subunit beta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rat
Immunogen Description	21-132AA
Other Names	Thyroid-stimulating hormone subunit beta,Thyrotropin beta chain,TSH-B
Accession No.	P04652Gene name:Tshb
Uniprot	P04652
GeneID	25653;
Calculated MW	12.32
Tag Info	His
Formulation	50mM NaH <sub>2</sub> PO <sub>4</sub> , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

## Background

Indispensable for the control of thyroid structure and metabolism.

## References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only