

Myoglobin

Catalog No: #AP78534



Package Size: #AP78534-1 50ug #AP78534-2 100ug #AP78534-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Myoglobin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	2-154AA
Other Names	null
Accession No.	P02144Gene name:MB
Uniprot	P02144
GeneID	4151;
Calculated MW	16.83
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole, 10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Serves as a reserve supply of oxygen and facilitates the movement of oxygen within muscles.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only