

## Endothelin-1

Catalog No: #AP78535


Package Size: #AP78535-1 50ug #AP78535-2 100ug #AP78535-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Endothelin-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	18aa-212aa
Other Names	Preproendothelin-1,PPET1
Accession No.	P05305Gene name:EDN1
Uniprot	P05305
GeneID	1906;
Calculated MW	21.45
Tag Info	His
Formulation	50mM NaH <sub>2</sub> PO <sub>4</sub> , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

## Images



## Background

Endothelins are endothelium-derived vasoconstrictor peptides.

## References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

---

Note: This product is for in vitro research use only