

NT-proBNP

Catalog No: #AP78536

Package Size: #AP78536-1 50ug #AP78536-2 100ug #AP78536-3 1mg


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	NT-proBNP
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	27-102aa
Other Names	N-terminal prohormone of brain natriuretic peptide
Accession No.	Gene name:NT-proBNP
Calculated MW	8.36
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.) Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Images



Background

The N-terminal prohormone of brain natriuretic peptide (NT-proBNP) is a 76 amino acid N-terminal fragment of brain natriuretic peptide.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only