

Hemoglobin subunit alpha

Catalog No: #AP78543



Package Size: #AP78543-1 50ug #AP78543-2 100ug #AP78543-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Hemoglobin subunit alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	Natural protein
Other Names	Alpha-globin,Hemoglobin alpha chain
Accession No.	P69905Gene name:HBA1
Uniprot	P69905
GeneID	3039;3040;
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Involved in oxygen transport from the lung to the various peripheral tissues.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only