

Phosphatidylinositol-glycan biosynthesis class F protein

Catalog No: #AP78544

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP78544-1 50ug #AP78544-2 100ug #AP78544-3 1mg

Description

Product Name	Phosphatidylinositol-glycan biosynthesis class F protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	1-180AA
Other Names	GPI11 homolog
Accession No.	Q07326Gene name:PIGF
Uniprot	Q07326
GeneID	5281;
Calculated MW	19.8
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.) Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Involved in GPI-anchor biosynthesis through the transfer of ethanolamine phosphate to the third mannose of GPI.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only