

## Neuritin

Catalog No: #AP78547

Package Size: #AP78547-1 50ug #AP78547-2 100ug #AP78547-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Neuritin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	28-116AA
Other Names	NRN
Accession No.	Q9NPD7Gene name:NRN1
Uniprot	Q9NPD7
GeneID	51299;
Calculated MW	9.79
Tag Info	His
Formulation	50mM NaH <sub>2</sub> PO <sub>4</sub> , 500mM NaCl Buffer with 500mM Imidazole, 10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

## Background

Promotes neurite outgrowth and especially branching of neuritic processes in primary hippocampal and cortical cells.

## References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only