

Neuritin

Catalog No: #AP78548



Package Size: #AP78548-1 50ug #AP78548-2 100ug #AP78548-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Neuritin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Mouse
Immunogen Description	28-116AA
Other Names	Cpg15,Candidate plasticity gene 15 protein
Accession No.	Q8CFV4Gene name:Nrn1
Uniprot	Q8CFV4
GeneID	68404;
Calculated MW	9.79
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Promotes neurite outgrowth and especially branching of neuritic processes in primary hippocampal and cortical cells.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only