

Fibroblast growth factor 2

Catalog No: #AP78555


Package Size: #AP78555-1 50ug #AP78555-2 100ug #AP78555-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Fibroblast growth factor 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Mouse
Immunogen Description	10-154aa
Other Names	Fgf-2,Basic fibroblast growth factor,Heparin-binding growth factor 2,bFGF,HBGF-2
Accession No.	P15655Gene name:Fgf2
Uniprot	P15655
GeneID	14173;
Calculated MW	15.95
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Images



Background

Plays an important role in the regulation of cell survival, cell division, angiogenesis, cell differentiation and cell migration. Functions as potent mitogen in vitro. Can induce angiogenesis.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only