

Collagen alpha-1(I) chain

Catalog No: #AP78568



Package Size: #AP78568-1 50ug #AP78568-2 100ug #AP78568-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Collagen alpha-1(I) chain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Mouse
Immunogen Description	23-151AA
Other Names	Cola1,Alpha-1 type I collagen
Accession No.	P11087Gene name:Col1a1
Uniprot	P11087
GeneID	12842;
Calculated MW	14.19
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Type I collagen is a member of group I collagen (fibrillar forming collagen).

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only