

N-terminal propeptide of Collagen alpha-1(III) chain

Catalog No: #AP78569



Package Size: #AP78569-1 50ug #AP78569-2 100ug #AP78569-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	N-terminal propeptide of Collagen alpha-1(III) chain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	103-206AA
Other Names	null
Accession No.	Gene name:PIIINP
Calculated MW	11.44
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.) Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only