

Somatostatin

Catalog No: #AP78582



Package Size: #AP78582-1 50ug #AP78582-2 100ug #AP78582-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Somatostatin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	AGCKNFFWKTFTSC
Other Names	Growth hormone release-inhibiting factor
Accession No.	P61278Gene name:SST
Uniprot	P61278
GeneID	6750;
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Somatostatin inhibits the release of somatotropin.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only