

Somatostatin

Catalog No: #AP78583

Package Size: #AP78583-1 50ug #AP78583-2 100ug #AP78583-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Somatostatin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Mouse
Immunogen Description	35-88AA
Other Names	Smst
Accession No.	P60041Gene name:Sst
Uniprot	P60041
GeneID	20604;
Calculated MW	5.94
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Somatostatin inhibits the release of somatotropin.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only