

Disintegrin and metalloproteinase domain-containing protein 8

Catalog No: #AP78596

Orders: order@signalwayantibody.com

Package Size: #AP78596-1 50ug #AP78596-2 100ug #AP78596-3 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Disintegrin and metalloproteinase domain-containing protein 8
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	677-824AA
Other Names	MS2,Cell surface antigen MS2,CD156a
Accession No.	P78325Gene name:ADAM8
Uniprot	P78325
GeneID	101;
Calculated MW	16.28
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Possible involvement in extravasation of leukocytes.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only