

Mucin-5B

Catalog No: #AP78615

Package Size: #AP78615-1 50ug #AP78615-2 100ug #AP78615-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mucin-5B
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	424-638AA
Other Names	MUC5,Cervical mucin,High molecular weight salivary mucin MG1,Mucin-5 subtype B, tracheobronchial,Sublingual gland mucin
Accession No.	Q9HC84Gene name:MUC5B
Uniprot	Q9HC84
GenelD	727897;
Calculated MW	23.65
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Gel-forming mucin that is thought to contribute to the lubricating and viscoelastic properties of whole saliva and cervical mucus.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only