

## C-C motif chemokine 24

Catalog No: #AP78622

Package Size: #AP78622-1 50ug #AP78622-2 100ug #AP78622-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	C-C motif chemokine 24
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Mouse
Immunogen Description	24-119aa
Other Names	Scya24, Eosinophil chemotactic protein 2, Eotaxin-2, Small-inducible cytokine A24
Accession No.	Q9JKC0 Gene name: Ccl24
Uniprot	Q9JKC0
GeneID	56221;
Calculated MW	10.56
Tag Info	His
Formulation	50mM NaH <sub>2</sub> PO <sub>4</sub> , 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.) Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

## Background

Chemotactic for resting T-lymphocytes, and eosinophils. Has lower chemotactic activity for neutrophils but none for monocytes and activated lymphocytes. Binds to CCR3.

## References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only