

# Thioredoxin

Catalog No: #AP78626



Package Size: #AP78626-1 50ug #AP78626-2 100ug #AP78626-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Thioredoxin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	Recombinant protein
Other Names	TRDX,TRX,TRX1,ATL-derived factor,Surface-associated sulphhydryl protein,ADF,SASP
Accession No.	P10599Gene name:TXN
Uniprot	P10599
GeneID	7295;
Tag Info	His
Formulation	50mM NaH <sub>2</sub> PO <sub>4</sub> , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

## Background

ADF augments the expression of the interleukin-2 receptor TAC (IL2R/P55).

## References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only