

Serpine A12

Catalog No: #AP78630



Package Size: #AP78630-1 50ug #AP78630-2 100ug #AP78630-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Serpine A12
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	21-414AA
Other Names	OL-64,Visceral adipose tissue-derived serine protease inhibitor,Visceral adipose-specific serpin,Vaspin
Accession No.	Q8IW75Gene name:SERPINA12
Uniprot	Q8IW75
GeneID	145264;
Calculated MW	43.34
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Adipokine that modulates insulin action by specifically inhibiting its target protease KLK7 in white adipose tissues.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only